

**Eye Controlled Speech  
Generating Device**  
**Your eyes say it all**



## MyTobii P10 Speech Generating Device



### Communication for millions

MyTobii P10 makes communication much more effective for users who previously relied on slow and cumbersome access methods such as single or multiple switch, scanning or head-mounted pointing devices.

### Who is helped by MyTobii?

User groups that benefit from MyTobii include people of all ages and cognitive functioning with Cerebral Palsy, ALS, Multiple Sclerosis, spinal cord injuries and many other conditions.

MyTobii's unique ability to deal with large head movements has opened up eye control for people who are unable to use other eye controlled systems, including those with CP.

### Easy to use

MyTobii takes only minutes to set up. The user does not need to "do" or "wear" anything. You just sit in front of the unit and look at keys or symbols displayed on the screen, as in front of a regular computer.

## The MyTobii P10 Speech Generating Device enables individuals with special needs to communicate using only their eyes.

After a fast, one-time calibration, no additional adjustments are needed. Eye control is fully automatic and so easy that individuals of all cognitive levels can communicate with the system.

### Reliable technology

The MyTobii P10 can be relied upon to perform. It operates with a high level of precision in nearly all light conditions. All different users, regardless of ethnic origin, can be tracked – including those with glasses or contact lenses. You can move freely in front of the screen without having to worry about interrupted tracking or reduced functionality.

### Designed for real life

The 15" screen and eye control functions are integrated into one sturdy unit. There are no fragile external cameras or lighting units.

It can be mounted for use at a desk, wheelchair, bed or anywhere suitable for the user. Just connect to a power source such as a wall socket, power wheelchair or separate battery. Robust casings will survive real-world use.

### Individual solutions for individual needs

MyTobii P10 is easily adjustable. The interface can be configured for many different skill levels to meet the needs of different user groups. Try-out tools provide help and guarantee optimal settings for each individual.

### Communication software of choice

MyTobii P10 includes Tobii Communicator Standard edition package, which allows for communication via text, or with over 10,000

SymbolStix® symbols. Tobii Communicator Premium, an upgrade, includes email, text messaging, and environmental control.

### Approved by Medicare

MyTobii P10 Basic is approved by Medicare as a dedicated SGD. Please contact Tobii or our local representative for information on how funding can be obtained.

### Product specifications MyTobii P10

Display	15" touch screen
Dimensions (WxHxD)	14.2" x 14.4" x 4.2"
Weight	11.7 lbs
Speakers	Integrated
Working distance	20-28 in
Freedom of head movement	12 x 6 x 8 in
Top head-motion speed	6 in/sec
Gaze data rate	40 Hz
Accuracy	0.5° (~0.2 in)
Max head movement compensation data	< 1 degree
Max long term deterioration, compensation error	< 1 degree
Computer	Integrated 1.5 GHz, Intel, 1 GB RAM, 80 GB Hard drive
Language	US English
Interface	2x USB II, 1x 100 MB LAN, 2x Joy switch connections
Mounting systems supported	Optional brackets for Vesa, Daessy and Rehadapt

### Accessories

The MyTobii P10 comes bundled with all that is needed to get started. A complete list of accessories is available at [www.tobiiATI.com](http://www.tobiiATI.com).

© Tobii®. Illustrations and specifications do not necessarily apply to products and services offered in each local market. Technical specifications are subject to change without prior notice. All other trademarks are the property of their respective owners.

#### NORTH AMERICA

Tobii ATI  
333 Elm Street  
Dedham, MA 02026  
Toll-Free: 800-793-9227  
Phone: 781-461-8200  
Fax: 781-461-8213  
Web: [www.tobiiati.com](http://www.tobiiati.com)  
[sales.us@tobii.com](mailto:sales.us@tobii.com)

#### EMEA

Tobii Technology AB  
Karlsrovägen 2D  
S-182 53 Danderyd  
Sweden  
Phone: +46 8 663 69 90  
Fax: +46 8 30 14 00  
[sales@tobii.com](mailto:sales@tobii.com)

#### ASIA

Tobii Technology, Ltd.  
3-4-13 Takanawa, Minato-ku  
Tokyo 108-0074  
Japan  
Phone: +81-3-5793-3316  
Fax: +81-3-5793-3317  
[sales.jp@tobii.com](mailto:sales.jp@tobii.com)

#### CENTRAL EUROPE

Tobii Technology GmbH  
Niedenau 45  
D-60325 Frankfurt am Main  
Germany  
Phone: +49 69 24 75 03 40  
Fax: +49 69 24 75 03 429  
[sales.de@tobii.com](mailto:sales.de@tobii.com)